



TECHNICAL SHEET

Product designation	MEDIBASE BLUE PLASTIC-COATED PAPER EXAMINATION SHEETS (6)
Reference	3173179
EAN	
Product image	
Product description	Blue plastic-coated paper examination sheets, for consultation tables. Highly resistant and waterproof, it has two material sides for optimum comfort of use.
Characteristics and compositions	<ul style="list-style-type: none">▪ 1 ply blue wadding paper laminated with polyethylene, embossed▪ Assembled with 4 knurls▪ Weight: 38 g/m² (28 g/m² wadding + 10 g/m² PE)▪ Width: 50 cm▪ Length: 57 m or 150 formats of 38 x 50 cm (±5%)▪ Pre-cut: 38 cm▪ Core inside diameter: 39 mm▪ Roll outside diameter: 9.5 < x < 9.8 cm
Packaging	Box of 6 rolls with carrying handles (gross weight: 6.5 kg)
Lifetime	/
Indications for use	
Legal information	Non-class product Made in France

Version 1.1
Update 2024/06/03 by CD