



Designation	POWDER-FREE LATEX EXAMINATION GLOVES ALOE VERA																							
Reference	XS – 5/6	S – 6/7	M – 7/8	L – 8/9																				
	3166225	3166226	3166227	3166228																				
Visual product																								
Product description	Powder-free latex examination gloves. Non-sterile, coated with Aloe Vera, ambidextrous, rolled cuff.																							
Characteristics and composition	<table border="1"> <thead> <tr> <th></th> <th>XS-5/6</th> <th>S-6/7</th> <th>M-7/8</th> <th>L-8/9</th> </tr> </thead> <tbody> <tr> <td>Width of the palm (mm) ± 10</td> <td>75</td> <td>80</td> <td>95</td> <td>111</td> </tr> <tr> <td>Single wall thickness (finger)</td> <td colspan="4">0.08 mm</td> </tr> <tr> <td>Single wall thickness (palm)</td> <td colspan="4">0.08 mm</td> </tr> </tbody> </table>					XS-5/6	S-6/7	M-7/8	L-8/9	Width of the palm (mm) ± 10	75	80	95	111	Single wall thickness (finger)	0.08 mm				Single wall thickness (palm)	0.08 mm			
		XS-5/6	S-6/7	M-7/8	L-8/9																			
Width of the palm (mm) ± 10	75	80	95	111																				
Single wall thickness (finger)	0.08 mm																							
Single wall thickness (palm)	0.08 mm																							
	<ul style="list-style-type: none"> Protein content: max 50 µg/dm². Powder level: max 2 mg/glove In compliance with standards: EN 455-1/2/3/4 																							
Packaging	Box included 100 gloves (± 2 gloves). 10 boxes per carton.																							
Shelf life	5 years from the date of manufacture.																							
Indication for use																								
Legal information	Medical Device of class I under MDR (EU) 2017/745 Manufacturer: MATANG MANUFACTURING SDN BHD.																							

Version 1.0
Update 2024/06/03 by CD