



Designation	Tray lining papers Medibase									
	Color	White	Blue	Green	Yellow	Orange	Pink	Lavender	Lime	Dark Blue
Reference(s)	Ref	116572	116744	116745	124643	3023338	3021948	3123064	3123065	3166481
Product picture										
Product description	Vapor and gas, permeable protections, which are suitable for sterilization and storage of instruments.									
Characteristics and composition	<ul style="list-style-type: none"> <li>- Dimensions: 18 x 28 cm</li> <li>- Very tear-resistant paper</li> <li>- High absorption capacity</li> <li>- Paper 55 g / m<sup>2</sup> (100% pure cellulose)</li> </ul>									
Packaging	Pack of 250 sheets, under plastic film protection									
Shelf life	N/A									
Indications of use										
Legal information	n.a									

Version 3

Updated 29/05/2024 by CD