



Product	Rubber cups	
Reference	Designation	Reference
	Bag of 144 cups	3135839
Picture		
Description	For polishing teeth after scaling.	
Characteristics & Technical data	<p>To use with a contra-angle. Have 6 intern lamellas: more contact for a better brightness. No dripping or splashes. Made of soft rubber to feet exactly to the tooth surface.</p> <p>10,000 max rpm</p>	
Packaging	Economic bag of 144 cups.	
Shelf life	5 years from the date of manufacture.	
Indication of use		
Legal informations	<p>Class IIA medical device. CE 1639 marking -SGS BELGIUM NV Notified body. Manufacturer: STODDARD & CO</p>	

Version 2 - 03/2025