



Designation	POWDER FREE NITRILE GLOVES															
<b>Reference(s)</b>	Color	XS	S	M	L	XL										
	Blue (200)	3135177	3135178	3135179	3135180	-										
	Blue (100)	3135172	3135173	3135174	3135175	3135176										
	White (100)	3135181	3135182	3135183	3135184	-										
	Lavender (100)	3135185	3135186	3135187	3135188	-										
<b>Visual</b>																
<b>Description</b>	Nitrile examination gloves, powder-free. Ambidextrous, finger-textured. Non sterile, beaded cuff															
<b>Characteristics and technical datas</b>	<ul style="list-style-type: none"> <li>• Excellent grip.</li> <li>• Allowing the manipulation of silicone.</li> <li>• Conforms to EN 455: 1-4 &amp; EN ISO 374-1;4;5.</li> <li>• Food contact approved</li> </ul>															
	<table border="1"> <thead> <tr> <th>Size</th> <th>Glove weight (g)</th> <th>Length (mm)</th> </tr> </thead> <tbody> <tr> <td><b>XS</b></td> <td>2.5 ± 0.2</td> <td rowspan="5">Min. 240</td> </tr> <tr> <td><b>S</b></td> <td>2.7 ± 0.2</td> </tr> <tr> <td><b>M</b></td> <td>2.9 ± 0.2</td> </tr> <tr> <td><b>L</b></td> <td>3.1 ± 0.2</td> </tr> <tr> <td><b>XL</b></td> <td>3.4 ± 0.2</td> </tr> </tbody> </table>		Size	Glove weight (g)	Length (mm)	<b>XS</b>	2.5 ± 0.2	Min. 240	<b>S</b>	2.7 ± 0.2	<b>M</b>	2.9 ± 0.2	<b>L</b>	3.1 ± 0.2	<b>XL</b>	3.4 ± 0.2
	Size	Glove weight (g)	Length (mm)													
	<b>XS</b>	2.5 ± 0.2	Min. 240													
	<b>S</b>	2.7 ± 0.2														
<b>M</b>	2.9 ± 0.2															
<b>L</b>	3.1 ± 0.2															
<b>XL</b>	3.4 ± 0.2															
<table border="1"> <thead> <tr> <th rowspan="3">Thickness (mm)</th> <th>Finger</th> <th>0.09</th> </tr> </thead> <tbody> <tr> <th>Palm</th> <th>0.05</th> </tr> <tr> <th>Cuff</th> <th>0.05</th> </tr> </tbody> </table>		Thickness (mm)	Finger	0.09	Palm	0.05	Cuff	0.05								
Thickness (mm)	Finger		0.09													
	Palm		0.05													
	Cuff	0.05														
<table border="1"> <thead> <tr> <th></th> <th>Before Aging</th> <th>After aging</th> </tr> </thead> <tbody> <tr> <td><b>Tensile strength (MPa)</b></td> <td>Min. 14</td> <td>Min. 14</td> </tr> <tr> <td><b>Elongation (%)</b></td> <td>Min. 500</td> <td>Min. 400</td> </tr> <tr> <td><b>Force at break (N)</b></td> <td>Min. 6</td> <td>Min. 6</td> </tr> </tbody> </table>			Before Aging	After aging	<b>Tensile strength (MPa)</b>	Min. 14	Min. 14	<b>Elongation (%)</b>	Min. 500	Min. 400	<b>Force at break (N)</b>	Min. 6	Min. 6			
	Before Aging	After aging														
<b>Tensile strength (MPa)</b>	Min. 14	Min. 14														
<b>Elongation (%)</b>	Min. 500	Min. 400														
<b>Force at break (N)</b>	Min. 6	Min. 6														
<table border="1"> <tbody> <tr> <td><b>Protein Level</b></td> <td>Not applicable</td> </tr> <tr> <td><b>Powder Level</b></td> <td>Max.2 mg/glove</td> </tr> </tbody> </table>		<b>Protein Level</b>	Not applicable	<b>Powder Level</b>	Max.2 mg/glove											
<b>Protein Level</b>	Not applicable															
<b>Powder Level</b>	Max.2 mg/glove															
<b>Packaging</b>	Box of 100 gloves Carton of 10 boxes															

<b>Storage</b>	Store protected from the light, water, humidity, and below 40°C
<b>Shelf life</b>	5 years since the date of manufacture
<b>Indication for use</b>	<p>EN ISO 374-5:2016    EN ISO 374-1:2016 Type B</p> 
<b>Legal information</b>	<p>Medical Device of class I under MDR REGULATION (EU) 2017/745</p> <p>PPE Cat III UNDER PPE REGULATION (EU) 2016/425 - CE Marking 0598</p>

Version 2.0 created by SMA

Date : 12/01/2023